

LIQUID CONTAINERS/PRESERVATIVES/HOLDING TIMES

All analysis needs to be kept cool to 4C except metals.

METAL ANALYSIS	CONTAINER	MAXIMUM HOLDING TIMES
Total ICP or ICP/MS Metals	HDPE preserved - 100ml	6 Months
Dissolved ICP or ICP/MS Metals(lab to filter)	HDPE unpreserved(lab to filter)-50ml	6 Months
Dissolved ICP or ICP/MS Metals(field filtered)	HDPE preserved(field filtered)-50ml	6 Months
Dissolved Mercury(lab to filter)	HDPE unpreserved(lab to filter)-50ml	28 Days
Dissolved Mercury(field filtered)	HDPE preserved HNO3(field filtered)-50ml	28 Days
Total Mercury	HDPE preserved HNO3 - 50ml	28 Days
Dissolved Low Level Mercury(field filtered)	125ml Amber unpreserved-50ml	28 Days
Dissolved Low Level Mercury(lab to filter)**24 Hours to filter	125ml Amber unpreserved-125ml	28 Days**
Total Low Level Mercury	125ml Amber unpreserved-125ml	28 Days
Hexavalent Chromium	HDPE unpreserved - 60ml	24 Hours
Ferrous Iron	HDPE unpreserved - 50ml	Immediate
ORGANIC ANALYSIS		
Alcohols	VOA unpreserved - 1 voa	7 Days
BNA	Amber unpreserved - 500ml	7 Days
DRO	Amber unpreserved - 500ml	7 Days
Formaldehyde	Amber unpreserved - 500ml	3 Days
Glycols	VOA unpreserved - 1 VOA	7 Days
GRO	VOA HCl - 2 VOA	7 Days
Herbicide	Amber unpreserved - 500ml	7 Days
PCB/Pesticide	Amber unpreserved - 500ml	7 Days
Phenols	Amber unpreserved - 500ml	7 Days
Phthalates	Amber unpreserved - 500ml	7 Days
PNA	Amber unpreserved - 500ml	7 Days
TPH	VOA HCl - 2 VOA	14 Days
Volatiles*1	VOA HCl - 2 VOA(preserved)	14 Days
MICROBIOLOGICAL		
Coliform Total+E.Coli(Drinking Water)	Pre-sterilized HDPE w/Sodium Thiosulfate	30 Hours
Coliform, Fecal	Pre-sterilized HDPE w/Sodium Thiosulfate	8 Hours
Coliform Total+E.Coli(Non Drinking Water)	Pre-sterilized HDPE w/Sodium Thiosulfate	8 Hours
Heterotropic Plate Count	Pre-sterilized HDPE w/Sodium Thiosulfate	8 Hours
INORGANIC ANALYSIS		
Alkalinity	HDPE unpreserved - 50ml	14 Days
Ammonia	HDPE H2SO4 - 10 ml	28 Days

BOD	HDPE unpreserved - 500ml	48 Hours
CBOD	HDPE unpreserved - 500ml	48 Hours
COD	HDPE H2SO4 - 10 ml	28 Days
Chloride	HDPE unpreserved - 20ml	28 Days
Chlorine(Residual)	HDPE unpreserved - 20ml	15 Minutes
Cyanide(Amenable)	HDPE NaOH - 150ml	14 Days
Cyanide(Available)*2	HDPE NaOH - 100ml	14 Days
Cyanide(Reactive)	Glass unpreserved - 50ml	ASAP
Cyanide(Total)	HDPE NaOH - 100ml	14 Days
Dissolved Oxygen	HDPE unpreserved - 500ml	15 Minutes
EOX/TX/TOX	Glass unpreserved - 100ml	28 Days
Flashpoint	Glass unpreserved - 100ml	ASAP
Fluoride	HDPE unpreserved - 20ml	28 Days
MBAS	Glass unpreserved - 300ml	48 Hours
Methane	VOA HCl - 2 VOA	14 Days
Nitrate	HDPE unpreserved - 20ml	48 Hours
Nitrite	HDPE unpreserved - 20ml	48 Hours
Oil & Grease	Glass H2SO4 - 1L	28 Days
pH(wastewater/drinking water)	HDPE unpreserved - 30ml	15 Minutes
pH(SW846)	HDPE unpreserved - 30ml	ASAP
Phenol(4AAP)	Glass H2SO4 - 400ml	28 Days
Phosphate (Ortho)	HDPE unpreserved - 10ml	48 Hours
Phosphate(Total)	HDPE H2SO4 - 25ml	28 Days
Phosphorus(Total)	HDPE H2SO4 - 25ml	28 Days
Resistivity	HDPE unpreserved - 30ml	28 Days
Specific Conductivity	HDPE unpreserved - 30ml	28 Days
Sulfates	HDPE unpreserved - 20ml	28 Days
Sulfide(Hydrogen Sulfide)*3	HDPE NaOH+ZnAc - 60ml	7 Days
Sulfide(Reactive)	Glass unpreserved - 50ml	ASAP
Sulfide(Total)	HDPE NaOH+ZnAc - 50ml	7 Days
Sulfite	HDPE unpreserved - 100ml	15 Minutes
TDS	HDPE unpreserved - 150ml	7 Days
TSS	HDPE unpreserved - 150ml	7 Days
TKN	HDPE H2SO4 - 50ml	28 Days

TOC	VOA HCl - 1 VOA	28 Days
Turbidity	HDPE unpreserved - 50ml	48 Hours

Footnotes

- *1. All volatile analysis require zero headspace.
- *2. Available Cyanide needs to be field filtered. Call lab for kit.
- *3. Hydrogen Sulfide - requires 2 separate containers(unpreserved HDPE&HDPE NaOH+ZnAc)

MS/MSD for any analysis requires more sample.

Abbreviations:

DW=Drinking Water
G=Glass
H2SO4=Sulfuric Acid
HCl=Hydrochloric Acid

HDPE=High Density Polyethylene
HNO3=Nitric Acid
NaOH=Sodium Hydroxide
ZnAc=Zinc Acetate